

**NOWdiagnostics, Inc.**<sup>™</sup>

1200 Stewart Place  
Springdale, AR 72764  
Phone: 479.966.4530  
Fax: 479.966.4631

**nowdx.com**

8001-PI  
Version 2.0

CE



**ADEXUSDX**<sup>™</sup>

**hCG TEST**

Product Family No. 8001 | For *in vitro* diagnostic use only

## TABLE OF CONTENTS

---

Intended Use .....	2
Features & Benefits .....	3
Summary & Explanation .....	4
Principle of the Procedure .....	4
Storage & Stability .....	6
Directions for Use .....	6
Materials .....	7
Sample Collection & Handling .....	9
Test Method .....	10
Interpretation of Results .....	14
Quality Control .....	16
Expected Results .....	17

Limitations .....	18
Disposal .....	19
Performance Characteristics .....	19
References .....	29

## INTENDED USE

---

**The ADEXUS<sup>Dx</sup>™ hCG Test** is an immunoassay used for the qualitative detection of human chorionic gonadotropin in human whole blood, plasma, or serum and is indicated as an aid for healthcare professionals in the diagnosis of early pregnancy.

## FEATURES

---

- Clinical relevant cutoff: 10 mIU/mL
- Small sample volume test (about 40  $\mu$ L)
- Simple one-step device
- Human capillary blood, plasma, or serum test

## BENEFITS

---

- Quickly analyze patients after 10 minutes
- Easy to use
- Accurate results (correlates with laboratory method)

## SUMMARY & EXPLANATION

---

Pregnancy tests are based on the detection of human chorionic gonadotropin (hCG), a hormone produced by the placenta around the fourth day after conception. hCG levels rise rapidly and double approximately every two days. The **ADEXUS**Dx™ hCG Test can detect hCG in a single drop of whole blood, serum, or plasma without additional reagents or materials.

## PRINCIPLE OF THE PROCEDURE

---

The **ADEXUS**Dx™ hCG Test is a rapid chromatographic immunoassay. A sample is dispensed at the Sample Application Zone to fill the Receiving Channel. When enough sample is in the Receiving Channel, the sample flows into a dry porous test strip composed of a plasma-separating membrane and an analytical membrane. The sample first passes through

## PRINCIPLE OF THE PROCEDURE (CONT.)

---

the plasma-separating membrane containing detector monoclonal mouse anti-hCG antibodies labeled with color. The detector antibodies bind to hCG in the sample to form soluble colored hCG complexes which then move to the analytical membrane. The colored hCG complexes are captured by polyclonal goat anti-hCG antibodies that are immobilized in the test band region of the analytical membrane. The appearance of a visible band indicates the sample contains a detectable level of hCG.

Excess detector antibodies will flow past the test band region and bind to immobilized polyclonal rabbit anti-mouse antibodies in the control band region of the analytical membrane. A visible control band with an absent test band assures that the negative result was not due to improper procedure.

## STORAGE & STABILITY

---

The **ADEXUS**Dx™ hCG Test is stable unopened at 4°C-30°C (39°F-86°F) until the expiration date.

## DIRECTIONS FOR USE

---

### Warnings and Precautionary Statements:

- For *in vitro* diagnostic use only.
- Dispose of expired tests and completed tests as biohazardous material.
- Do not use expired tests.
- Do not use if the pouch seal is compromised.
- Single use only.

## DIRECTIONS FOR USE (CONT.)

---

**CAUTION:** The device contains material of human or animal origin and should be handled as a potential carrier and transmitter of disease.

## MATERIALS

---

### MATERIAL PROVIDED:

- The **ADEXUS**Dx™ hCG Test is packaged in individually sealed pouches.

### MATERIAL REQUIRED BUT NOT PROVIDED:

#### Capillary Blood Test

- Commercially available lancet for obtaining blood from a finger puncture
- Alcohol prep pads
- Timer

#### Venous Blood, Plasma, and Serum Test

- Heparinized sample collection container for venous blood
- Appropriate sample collection container and equipment for obtaining plasma or serum
- Sample transfer device
- Timer

## SAMPLE COLLECTION & HANDLING

### Capillary Blood Test

- A blood sample from a finger puncture should be obtained immediately before application to the test.



### Venous Blood, Plasma, and Serum Test

- For venous blood or plasma tests, obtain blood by venipuncture. In a heparinized container, collect blood for the test or process the blood to obtain plasma.
- For serum tests, obtain blood by venipuncture. In an appropriate container, collect blood and process the blood to obtain serum.

Venous blood samples may be stored at room temperature 15°C-25°C, (59°F-77°F) for up to 4 hours prior to testing. Plasma and serum samples may be stored for 72 hours at 2°C - 8°C (36°F - 46°F) prior to testing.

9 Samples should be at room temperature for the test.

## TEST METHOD



## TEST METHOD (CONT.)

---

### For Capillary Whole Blood Samples

1. Remove the **ADEXUS**Dx™ hCG Test from the foil pouch and discard the desiccant and packaging.
2. To promote blood flow to the fingertips, ask the person to lower the hand to be punctured and massage the fingers toward the tips of the fingers.
3. Choose either the “ring” or “middle” finger for the test site. Wipe the fingertip with an alcohol prep pad and allow the finger to thoroughly air-dry.
4. Prepare the lancet following the manufacturer’s instructions.
5. Select the fingertip test site. Do not use the center or top of the finger, because these are the most sensitive areas of the finger.
6. Squeeze the center of the finger (with one hand) at the same time as pressing firmly on the test site with the lancet (with the other hand).
7. Release the lancet following the manufacturer’s instructions.
8. Gently massage the fingertip toward the puncture site to obtain a full drop of blood.
9. Touch the puncture site with the Sample Application Zone of the **ADEXUS**Dx™ hCG Test.
10. Wait until all blood sample is drawn into the **ADEXUS**Dx™ hCG Test. If the blood sample has not reached the Fill Line, continue to massage the finger to obtain more blood until the blood sample reaches the Fill Line.
11. Start the timer.
12. After 10 minutes, observe the appearance of bands in the Test Window at the Control Position **(C)** and the Test Position **(T)**. See Interpretation of Results. Do not interpret the results of the test after 30 minutes.


## TEST METHOD (CONT.)

### For Venous Blood, Plasma, or Serum Samples

1. Remove the **ADEXUS**Dx™ hCG Test from the foil pouch and discard the desiccant and packaging.
2. Apply the sample to the **ADEXUS**Dx™ hCG Test Sample Application Zone using a liquid transfer device (not supplied) until the sample reaches the Fill Line (about 40µL).
3. Start the timer.
4. After 10 minutes, observe the appearance of bands in the Test Window at the Control Position (**C**) and the Test Position (**T**). See Interpretation of Results. Do not interpret the results of the test after 30 minutes.

## INTERPRETATION OF RESULTS

A band of any signal intensity appearing at the Test Position (T) is considered positive for hCG.



## INTERPRETATION OF RESULTS (CONT.)

---

### Notes:

- The control band and the test band may vary in color intensity.
- The intensity of the test band will vary depending on the concentration of hCG present in the specimen. However, the quantitative value of hCG cannot be determined by this qualitative test.
- The color intensity of the bands will increase slowly with time as a result of sample evaporation. Therefore, do not read the test after 30 minutes.

## QUALITY CONTROL

---

The **ADEXUS**Dx™ hCG Test has an internal quality control band. In cases where the test band is absent, the appearance of the control band assures that the sample was applied correctly and that proper chromatography in the test occurred.

It is recommended that external quality control materials be used to verify the performance of the **ADEXUS**Dx™ hCG Test. A negative hCG control (< 5 mIU hCG/mL) and a positive hCG control (> 10 mIU hCG/mL) should be used. It is suggested that federal, state, and local guidelines be followed.

## EXPECTED RESULTS

---

hCG in nonpregnant women is normally not detected by the **ADEXUS**Dx™ hCG Test. hCG concentrations are generally between 5 and 50 mIU/mL within 1 week of gestational age or 3 weeks after the last menstrual period (LMP). hCG levels rise rapidly, doubling about every two days, and peak to greater than 100,000 mIU/mL in the latter part of the first trimester of pregnancy. If a test result is negative and pregnancy is suspected, repeat the test after 2 days or more.

A faintly positive test can be confirmed by repeating the test after 2 days or more.

## LIMITATIONS

---

- Other clinical conditions that produce hCG include trophoblastic diseases and germ cell tumors.
- Highly sensitive pregnancy tests such as the **ADEXUS**Dx™ hCG Test may detect pregnancies ending in an early pregnancy loss, which occurs in about one out of every four pregnancies. Results should be used in conjunction with other clinical and laboratory data.
- Very strong positive samples may cause the control band to become faint. A high dose hook effect was not evident at 150,000 mIU hCG/mL.
- As with any other assay using mouse antibodies, the possibility exists for interference by human anti-mouse antibodies (HAMA) or other interfering substances in the sample.

## DISPOSAL

---

Dispose of the expired tests and completed tests as biohazardous material.

## PERFORMANCE CHARACTERISTICS

---

### Sensitivity:

The test band was designed to be visible when a specimen containing an hCG concentration of approximately 10 mIU/mL is analyzed by the test. The test was calibrated with the Abbott Architect for total  $\beta$ -hCG which is calibrated with material referenced to WHO, 3<sup>rd</sup> International Standard 75/537.

### Assay Precision Near the Cutoff:

Three hCG positive patient serum pools were prepared to have concentrations near 20, 10, and 5 mIU hCG/mL. Twenty (20) replicates of each pool for a total of 60 samples were tested with the **ADEXUS**Dx™ hCG Test. The **ADEXUS**Dx™ hCG Test identified the 40 positive samples (20 and 10 mIU hCG/mL) and the 20 negative samples (5 mIU hCG/mL).

### Method Comparison:

A comparison study between the **ADEXUS**Dx™ hCG Test and another commercially available serum hCG assay was performed. Fifty-nine (59) and 57 female serum samples with hCG concentrations above and below 10 mIU/mL, respectively, were tested.

		ADEXUSDx™ hCG Test	
		Positive	Negative
Other Assay	Positive	59	4
	Negative	0	53

For hCG-positive serum samples (10-200 mIU/mL), there was 100% agreement between the two assays. For normal serum samples, there were four discrepancies in which the other assay gave positive results and the **ADEXUSDx™** hCG Test gave negative results. The discrepant samples had hCG concentrations in the range of 2-6 mIU/mL, and were correctly reported as negative by the **ADEXUSDx™** hCG Test based on the **ADEXUSDx™** hCG Test 10 mIU hCG/mL cutoff value. There was 100% agreement of the **ADEXUSDx™** hCG Test with the other assay for the remaining normal samples.

### Proficiency Testing:

Three samples of normal human whole blood were spiked with hCG to attain concentrations less than or above the cutoff value of 10 mIU/mL. Three different hospitals and on-site for a total of 4 sites were provided with 5 replicates of each sample. The replicates were sent as blind samples with a unique label for each of the 15 tubes. There was complete agreement of the results between sites.

Site	Number of Positive Results/Total			Agreement within Run
	Below cutoff (not spiked)	3x cutoff	6x cutoff	
1	0/5	5/5	5/5	100%
2	0/5	5/5	5/5	100%
3	0/5	5/5	5/5	100%
On-site	0/5	5/5	5/5	100%
<b>Agreement Between Sites</b>	<b>100%</b>	<b>100%</b>	<b>100%</b>	

### Dilution Study:

Three serial two-fold dilutions of hCG-positive plasma were made. Five replicates of each sample were applied to the **ADEXUS**Dx™ hCG Test. Samples with hCG concentrations above 10 mIU/mL all showed positive results, and below 10 mIU/mL all showed negative results.

hCG conc. (mIU/mL)	Number of Positive/Total
24.8	5/5
12.4	5/5
6.2	0/5
3.1	0/5

## Correlation Between Whole Blood and Plasma

Heparinized whole blood from a non-pregnant female was spiked with hCG to attain samples with concentrations <1.2, 4, 9, and 19 mIU hCG/mL when plasma recovered from the spiked blood was quantified. Both blood and the corresponding plasma samples for a total of 4 sets of samples were tested with the **ADEXUS**Dx™ hCG Test. Each set was tested 5 times.

Results from the blood and plasma test were correlated.

		Whole Blood	
		Positive	Negative
Plasma	Positive	10	0
	Negative	2	8

The agreement between whole blood and plasma is 100% for hCG values greater than or equal to 9 mIU/mL. Two tests gave a weak positive signal at 4 mIU hCG/mL for whole blood.

## Hook Effect:

Normal human serum was spiked with hCG to attain a concentration of 150,000 mIU hCG/mL. Five (5) replicates were tested with the **ADEXUS**Dx™ hCG Test. In all cases the test band signal was very strong.

## Capillary Blood Study:

Capillary blood samples from 10 non-pregnant female donors were taken using commercially available lancets to puncture fingers. The samples were applied to the **ADEXUS**Dx™ hCG Test and examined after 10 minutes. For all the tests, the control band was present, the red blood cell front was below the test window, the plasma front reached the top end of the test strip, and hemolysis did not occur. The test band was absent in all the tests indicating negative results.

**Specificity:**

The following related hormones were added to human negative control serum.

hTSH	1000 $\mu$ IU/mL
hLH	500 mIU/mL
hFSH	1000 mIU/mL

Negative results were obtained at the given concentrations.

**Interfering substances:**

The following potentially interfering substances were added to serum samples containing 0 and 10 mIU hCG/mL.

Protein	14 g/dL
Hemoglobin	250 mg/dL
Triglycerides	2000 mg/dL
Bilirubin	15 mg/dL

None of the substances at the concentration tested interfered with the **ADEXUS**Dx™ hCG Test.

## REFERENCES

---

1. Ho HH, O'Connor JF, Nakajima ST, Tieu J, Overstreet JW, Lasley BL. Characterization of human chorionic gonadotropin in normal and abnormal pregnancies. *Early Pregnancy*. 1997; 3(3):213-24.
2. Lenton EA, Neal LM, Sulaiman R. Plasma concentration of human chorionic gonadotropin from the time of implantation until the second week of pregnancy. *Fertil. Steril*. 1982; 37(6):773-778.
3. Batzer FR. Hormonal evaluation of early pregnancy. *Fertil. Steril*. 1980; 34(1):1-13.
4. Braunstein GD, Rasor J, Danzer H, Adler D, Wade ME. Serum human chorionic gonadotropin levels throughout normal pregnancy. *Am. J. Obstet. Gynecol*. 1976; 126(6):678-681.



Simple test



Perform in-office  
No outside lab required



Fast and  
accurate results

**ADEXUS**Dx™ is reinventing how diagnostics are performed. Our patented approach allows physicians to easily and accurately receive results quickly by using the tiniest blood sample. The focus of our rapid format for the laboratory and point-of-care diagnostics is to meet the needs of today's healthcare professionals and improve patient care.