

# UNAIDED AND AIDED LOCALIZATION PERFORMANCE

Denise M. Keenan, M.A. // Jenn Schumacher, B.A.  
Widex Office of Research in Clinical Amplification (ORCA)



## INTRODUCTION

The pinna attenuates sounds that originate from the rear by approximately 5 dB in the 2000-6000 Hz range. By positioning the microphone of a hearing aid outside the concha area, such as a behind the ear [BTE] instrument, this directional advantage from the pinna is lost. A new algorithm was developed to mimic the directional advantage of the pinna for BTE instruments [the Digital Pinna]. This algorithm would assist with front/back localization. This case study poster will examine the front/back localization abilities of four individuals tested in the unaided and aided conditions with and without the new algorithm.

## METHODS

### Hearing instrument

- Widex Clear 440 Passion, Receiver in Canal
  - 15 channel
  - Fully adaptive directional microphone
  - Speech intelligibility based noise reduction
  - Multi directional active feedback cancellation
  - Inter ear communication [compression, noise reduction, volume, program, zen tones...]
- Custom CAMISHA shells made for each participant
- Venting based on hearing loss at 500 Hz
- <20 dB = open fit; 20-29 dB = 1.5-2mm; 30-39 dB = 1-1.5mm; 40-49 dB = 0.5-1mm; 50-60 dB = 0-0.5mm
- For testing purposes
  - Digital Pinna feature turned on and off for test comparison
  - Hear aid set in omni directional mode
  - Noise reduction off
  - Feedback cancellation set to default
  - Inter ear communication off

### Set up

- Participants were seated in a 10'x10'x6'6" double wall sound treated booth.
- Participants faced a 0° speaker in a 12 speaker array [Figure 1].
  - Each speaker separated by 30°
  - Each speaker one meter from center
  - A 17" Planar PT 1700 MU touch screen was placed in front of listener for responses.
  - A foam headrest was used to remind the participant to keep head position straight.
  - A 300 ms high pass filtered noise [above 2000 Hz] was presented at 30 dB sensation level [SL] in sound field.
  - The test signal was presented randomly from one of the twelve speakers. The test ended when the signal had been presented a total of three times from each of the twelve speakers [36 total test signals].




Figure 1: Representation of speaker array

## METHODS [cont.]

### Procedures

- Visit 1:**
- Participant seen for unaided testing and insert earphone puretone audiometric check.
  - Ear impressions were taken and custom CAMISHA earmolds ordered.
- Visit 2:**
- Instruments programmed with in-situ threshold measurements [Sensogram] and feedback test.
  - Default gain settings were used. Fine tuning permitted if participant reported issue with loudness.
  - Localization testing conducted.
  - For home use, hearing aid settings included: Locator [automatic adaptive directional microphone] with Digital Pinna algorithm, Speech Enhancer [SII based noise reduction], feedback cancellation set to default, and inter-ear communication was active [noise reduction, compression, feedback cancellation ...].
- Visit 3:**
- One month after fitting. Localization testing repeated.
- Visit 4:**
- Two months after hearing aid fitting. Localization testing repeated.

## RESULTS

Even though the participants were evaluated at 12 different angles for each test, for reporting purposes of front and back localization, the front score was comprised of 330°, 0°, and 30°; and the back score was comprised of 150°, 180°, and 210°.

### Participant 1

- Participant 1: a 73 year old female
- Presented with a mild to moderate sensorineural hearing loss [Figure 2].
  - She had 2 years of hearing aid experience with a BTE instrument.
  - Her own hearing aid had automatic adaptive directional microphones.
- Results shown in Figure 3:**
- Unaided testing indicated that she had better localization to the front than to the back [by 44%].
  - With no Digital Pinna algorithm, she was unable to identify sounds to the back [which did not improve over time].
  - With the Digital Pinna algorithm, she was able to identify sounds to the front and sounds to the back equally [89% correct].
  - The results of this participant with the Digital Pinna were very similar to results obtained for those with normal hearing on this test.




Figure 2: Audiogram for Participant 1




Figure 3: Front/Back results for Participant 1

## RESULTS [cont.]

### Participant 2

- Participant 2: a 69 year old male
- Presented with a precipitous high frequency hearing loss [Figure 4].
  - He had 7 years of hearing aid experience with an ITE instrument.
  - His hearing aid only had an omni directional microphone.
- Results shown in Figure 5:**
- Unaided testing indicated that he had better localization to the front than to the back [by 44%].
  - With no Digital Pinna algorithm, he localized more to the back than to the front [by 55%]; this did not improve over time.
  - With the Digital Pinna algorithm, identification to the back [89%] as well as to the front [22%] improved over no Digital Pinna.
  - Front localization improved over time.
  - Localization ability to the back was considerably better with the Digital Pinna algorithm than unaided. After experience, localization to the front was similar to unaided.




Figure 4: Audiogram for Participant 2




Figure 5: Front/Back results for Participant 2

### Participant 3

- Participant 3: a 76 year old male
- Presented with a mild to moderate sensorineural hearing loss [Figure 6].
  - He had no hearing aid experience.
- Results shown in Figure 7:**
- Unaided testing indicated that he had better localization to the front than to the back.
  - With no Digital Pinna algorithm, he was unable to localize in either direction [which did not improve over time].
  - With the Digital Pinna algorithm, he began to localize some sounds to the back.
  - With experience, localization to the back improved considerably and he began to localize sounds in the front as well.
  - If given more time to acclimatize to amplification, it would be interesting to see if this participant would continue to improve localization ability.




Figure 6: Audiogram for Participant 3




Figure 7: Front/Back results for Participant 3

## RESULTS [cont.]

### Participant 4

- Participant 4: a 65 year old male.
- Presented with a precipitous high frequency hearing loss [Figure 8].
  - He had no hearing aid experience.
- Results shown in Figure 9:**
- Unaided testing indicated that he had minimal localization to the front and no localization to the back.
  - With no Digital Pinna algorithm, he was not able to localize sounds to the front and was able to localize some of the sounds to the back [33%]. This did not improve over time.
  - With the Digital Pinna algorithm, he was able to localize sounds to the back [67%] and began to localize sounds to the front [22%].
  - With experience, localization to the front started to improve and localization to the back decreased slightly by 10%.
  - As with Participant 3 with no hearing aid experience, it would be interesting to see if his localization skills continued to improve with time.




Figure 8: Audiogram for Participant 4




Figure 9: Front/Back results for Participant 4

## OBSERVATIONS

- Participants 1 and 3, although having similar hearing losses, had very different localization ability.
- Participants 2 and 4, also having similar hearing losses, had very different localization ability.
- The participants with hearing aid experience, Participants 1 and 2, had better localization ability with amplification than did those without experience, Participants 3 and 4.
- The participants without hearing aid experience did improve over time. It would be interesting to see if this improvement continued with more hearing aid experience.
- These participants had better localization ability with the use of amplification than their unaided localization.
- All four participants demonstrated better localization ability with the use of the Digital Pinna than without the use of the Digital Pinna algorithm.

## CONCLUSIONS

- Localization ability was not determined by auditory thresholds alone.
- Localization ability was influenced by hearing aid experience for these listeners.
- Further investigation will be needed to examine if localization ability continues to improve over time.
- Further investigation will be needed to examine why there are differences in localization ability for the hearing impaired besides auditory thresholds; such as damage to inner hair cells.