IKA

designed for scientists

IKA AOD 1.1

<i>Betriebsanleitung</i>	DE	2
<i>Operating instructions</i>	EN	10
<i>Mode d'emploi</i>	FR	18
<i>Indicaciones de seguridad</i>	ES	26
<i>Veiligheidsinstructies</i>	NL	28
<i>Norme di sicurezza</i>	IT	30
<i>Säkerhetsanvisningar</i>	SV	32
<i>Sikkerhedshenvisninger</i>	DA	34
<i>Sikkerhetsanvisninger</i>	NO	36
<i>Turvallisuusohjeet</i>	FI	38
<i>Υποδείξεις ασφάλειας</i>	EL	40
<i>Instruções de segurança</i>	PT	42

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
Explanation of symbols

This symbol identifies information **that is absolutely essential for your safety and health**. Failure to observe this information may result in negative effects for your health and injuries.

This symbol identifies information that is important for problem-free operation as well as for working with the decomposition vessel. Failure to observe these instructions may result in inaccurate measurement results or damage to the decomposition vessel.

Contents

Explanation of symbols	10
Contents	10
Safety Instructions	10
Included with delivery AOD 1.1	12
Opening the decomposition vessel	12
Fastening cotton thread	12
Combustion procedure	13
Cleaning the decomposition vessel	13
Halogens and decomposition	13
Care and maintenance	13
Decomposition vessel identification	15
Accessories	16
Technical data for AOD 1.1	17
Guarantee	17
Spare parts diagram for the AOD 1.1	17

Safety Instructions

Operation of decomposition vessels is only permitted in combination with the AOD 1 system.

Operating requirements:

The maximum amount of energy input into the decomposition vessel must not exceed 30000 J. (Select the weight of the sample accordingly). The permissible operating pressure of 195 bar must not be exceeded. The maximum permissible operating temperature must not exceed 50 °C.

Do not fill the decomposition vessel too full of the sample. Only fill the decomposition vessel with oxygen up to a maximum pressure of 40 bar. Monitor the adjusted pressure on the pressure reducer. Perform a check before every combustion to ensure there are no leaks (see page 13).

Explosive substances:

Many substances tend to combust in an explosive manner (for example because of the formation of peroxide). This may cause the decomposition vessel to burst.


The AOD 1.1 decomposition vessel must not be used for examinations on samples that are capable of exploding. It is absolutely essential to use a special high-pressure decomposition vessel to contain the sample in these cases!

Notes on the sample:

Substances of which the combustion behavior is not known must be examined for their combustion behavior before combustion in the decomposition vessel AOD 1.1 (danger of explosion). If you are burning unknown samples, leave the room or keep a safe distance between you and the calorimeter.

Benzoic acid must only be burned in the form of pellets! Combustible dust and powder must be compressed into pellets before combustion. Oven-dry dust and powder such as wood chips, hay, straw, etc. burn in an explosive manner! They must be moistened first! Readily combustible liquids with a low vapor pressure must not be come in direct contact with the cotton thread (for example tetramethyl dihydrogen disiloxan)!

Combustion residue, auxiliary materials

In addition, toxic residues of combustion are possible in the form of gasses, ash or precipitates on the inner wall of the decomposition vessel, for example.

Observe the accident prevention requirements applicable to the activity and the work station. Wear personal safety equipment.

When handling combustion samples, combustion residues and auxiliary materials, the appropriate safety requirements must be observed. The following are examples of substances that may cause dangers:

corrosive, easily flammable, capable of exploding
contaminated with bacteria, toxic

Oxygen

When working with oxygen, observe the appropriate requirements.

Danger warning: As a compressed gas, oxygen promotes combustion, supports combustion intensively and may react violently with combustible substances. **Do not use any oil or grease!**

Using a crucible made of stainless steel

When using crucibles made of stainless steel, their condition should be carefully checked after every experiment.

A reduction in the thickness of the material may cause the crucible to burn and may damage the AOD 1.1 decomposition vessel.

For reasons of safety, crucibles must not be used any more after a maximum of 25 combustion procedures.

Specification of the decomposition vessel:

The decomposition vessel is manufactured in accordance with the regulation for pressure vessels 97/ 23/ EC. The decomposition vessel has been subjected to a pressure test at a test pressure of 280 bar and a leak test with oxygen at 30 bar.

The AOD 1.1 decomposition vessel is an experiment autoclave and must be tested by a professionally trained person each time before they are used.

An individual application is understood here to mean a series of experiments that are performed under roughly the same conditions in terms of pressure and temperature Experiment autoclaves

must be operated in special chambers / protective device.

Repeated tests

The decomposition vessel must be subject to repeated tests (internal tests and pressure tests) by a **person with professional training**. The intervals between tests must be determined by the operator based on experience, operating manner and the material used in the decomposition vessel.

The guarantee loses its validity if mechanical modifications are made to the experiment autoclaves or if stability can no longer be guaranteed as a result of heavy corrosion (for example holes eaten in it by halogens).

The threading on the body of the decomposition vessel and cap screw in particular are subject to a high level of mechanical stress and must therefore be monitored regularly for wear and tear.

The condition of the seals must be checked for functionality must be ensured by means of a test for leaks.

(see the sections entitled "Care and maintenance; Leak test")

Pressure tests and service tasks on the decomposition vessel must only be performed by **persons with professional training**.

We recommend that the decomposition vessel be sent into our factory for inspection and repairs if necessary after either 1000 experiments or after one year or, depending on the application, even sooner than this.

Definition of person with professional training

A person with professional training as defined in these operating instructions is someone

1. whose training, knowledge and experience gained through practical activities ensures that that person will perform the tests in a proper manner.
2. who is sufficiently reliable
3. who is not subject to any instructions in terms of testing activity
4. who is equipped with suitable testing equipment if necessary
5. who can provide suitable proof demonstrating compliance with the requirements listed in 1.

Operating pressure containers

National regulations and laws for operating pressure containers must be observed!

Anyone who operates a pressure container must maintain it in proper condition, must monitor it and perform necessary maintenance and repair tasks without delay, and must take measures appropriate for the circumstances to ensure safety.

A pressure container must not be operated if it exhibits defects that could endanger those working with it or third parties.

You can obtain a copy of the pressure vessel regulation from Carl Heymann Press or Beuth Press.

Included with delivery of AOD 1.1


1 x decomposition vessel	1 x operating instructions
1 x ignition wire	1 x C 4 combustion crucible
1 x O-ring set	1 x SW5,5 double open ended wrench
1 x C710.4 cotton thread	1 x snap ring
1 x grip	

Opening the decomposition vessel

Unscrew union nut and use the handle in order to pull out the cover.


Unscrew union nut and use the handle in order to pull out the cover.

Fastening cotton thread

The decomposition vessel is equipped with a fixed ignition wire. A cotton thread has to be fastened to the ignition wire in order to carry out combustion.

Combustion procedure

For the decomposition vessel charging procedure with fuel samples and the general procedure for performing combustion please refer to the operating instructions for the IKA AOD 1 system.

Cleaning the decomposition vessel

After each combustion the condensate has to be cleaned from all of the parts of the decomposition vessel which come into contact with the combustion products. It suffices to wipe out the parts with an absorbent, non-fibrous cloth. Combustion residues in the crucible, e.g. carbon black or ash are removed in the same manner. The decomposition vessel may not be mechanically cleaned under any circumstances if it is not capable of being cleaned in the manner described (e.g. due to burn-in or pitting).

In such cases please contact your competent service center or send in the decomposition vessel to the factory for inspection and cleaning.

Halogens and decomposition

Decomposition vessel AOD 1.1 is obtainable in a modified form for determining halogens and sulfur. This means that the inner surfaces and the fittings of the decomposition vessel are catalytically activated by means of a patented procedure.

Parts thus treated are recognizable because of their opaque, tarnished and partly spotty surfaces. This appearance of this surface is intentional and should not be confused with traces of wear. This state is important for the following analysis and may not be treated with brushes or other hard materials.

Care and maintenance

The work described in the following may only be performed by personnel with professional training or IKA.

All work on the fittings of the decomposition vessel require a leak test and new pressure testing.

Furthermore, reference is made to the regulation for pressure vessels 97/23/EC.

Replacing the ignition wire

The AOD 1.1 decomposition vessel is equipped with a fixed ignition wire. Since the ignition wire is subjected to a high degree of material stress during combustion, it should be replaced after about 100 tests at the latest. However, this largely depends on the type of samples fired.

- Unscrew the upper hexagonal nut (item 11) on the ignition electrode (item 5) or the crucible holder (item 12) using a double open ended spanner SW 5,5 which is included in the delivery quantity.

Replacing the ignition electrode

The ignition electrode has to be replaced if it has become worn as a combustion in the decomposition vessel.

- Open the decomposition vessel as described under „Opening the decomposition vessel“.
- Take hold of the complete cover unit.
- Remove the ignition wire
- Unscrew the hexagonal nut (item 11) in a counterclockwise direction using the double open ended spanner (SW 5,5) while firmly holding the ignition electrode (item 5). Screw the hexagonal nut completely off of the ignition electrode and put it aside for later assembly.
- Pull the ignition electrode out of the cover. Make certain that the insulating bushing (item 8) and the bearing bush (item 9) are removed with the electrode. Clean both parts if required.

 **A new o-ring (item 7) should be used when the new electrode is inserted.**

Information on how to replace the O-ring is provided further on in this chapter. Make certain that there are no dirt particles in the cover when reassembling the unit. Use a paper towel to clean the sealing surfaces. Installation is performed in the reverse order.

- Put the rate insulating bushing (item 8 and 13) and the bearing bushing (item 9) on to ignition electrode.

- Insert

the new ignition electrode contained in the set of electrodes with deflector (item 2002) into the cover.


- The hexagonal nut (item 11) should only be slightly tightened with the spanner.

A leak test should be performed after assembly of the decomposition vessel has been completed.

- Fill the decomposition vessel with oxygen. Use the C 48 oxygen station in order to do so.

- Hold the decomposition vessel with the cover upside down in a vessel that is completely filled with water.

Be careful not to immerse only half of the decomposition vessel at most.

- Check the decomposition vessel for escaping bubbles of oxygen.

The leak test does not replace pressure testing !

Oxygen valve maintenance

The oxygen valve is used to fill and empty the decomposition vessel. It is a non-returnable valve that only closes when the decomposition vessel is full with oxygen.

The ignition electrode (item 5) serves at the same time as a valve piston. The sealing surface is in the cover of the decomposition vessel. The O-ring sealing wears through use and because of aggressive combustion gases, which means that the decomposition vessel may not be tight.

Check the sealings for wear in accordance with regular intervals stipulated by the safety provisions. If the decomposition leaks during combustion in the calorimeter the escaping oxygen may damage the decomposition vessel and the calorimeter.

- In order to replace the O-ring the ignition electrode must be removed from the cover as described.
- Remove the O-ring (item 7) by first pulling the bearing bushing (item 9) off of the ignition electrode and then removing the O-ring.
- Place the new O-ring from the set of O-rings onto the ignition electrode (item 5) from above and then push the bearing bushing (item 9) until it again reaches the ignition electrode.

Make certain that the O-ring does not become damaged by sharp edges.

After installation of the ignition electrode the bearing bushing is tightened with the hexagonal nut (item 11) in order to prevent it from slipping out of place.

Replacing the main sealing

The main sealing (item 4) is used to seal off the cover (item 2) from the bottom part (item 1).

The main sealing and the sealing surface on the bottom part must be checked for wear by an expert after each combustion operation. The main sealing must be replaced if wear is evident.

- Remove the main sealing (item 4) by using a suitable tool to pry the O-ring out of the groove in the cover (item 2).
- Pull the O-ring toward the side of the cover that is smaller in diameter.

Make certain not to damage the cover with tool. This may cause the decomposition vessel to leak. The sealing that has thus been replaced may not be reused. Replaced the old sealing with a new one from the set of O-rings.

- The new sealing is also mounted on the side of the cover that is smaller in diameter. Place the O-ring on one side the groove in the cover.
 - Use your fingers to pull the O-ring over the edge of the cover until it snaps into groove on its own.
- Perform a leak test after the decomposition vessel has been assembled.

- Fill the decomposition vessel with oxygen (C 48).
- Hold the decomposition vessel with the cover upside down in a vessel that is filled with water.
- Check the decomposition vessel for escaping bubbles of oxygen.

The leak test does not replace pressure testing !

Replacing the crucible holder

The crucible holder must be replaced if it has become worn as a result of the combustion operations.


- Open the decomposition vessel as described under „Opening the decomposition vessel“.
- Remove ignition wire and ignition electrode as described before.

- Unscrew the hexagonal nut adjacent to the inside cover of the crucible holder in a counterclockwise direction using the SW 5,5 double open ended spanner.
 - Screw the crucible holder out of the cover.
- Installation of the new crucible holder is performed in the reverse order. The crucible holder has to be aligned (see page 14) before tightening the hexagonal nut (item 14). **The hexagonal nut is tightened only slightly with the spanner.**

Decomposition vessel identification

The individual parts of the decomposition vessel are clearly identified by means of the engraved serial number and the technical data. The data is in accordance with the regulation for pressure device.

The individual parts of several decomposition vessels (cover, cap screw and bottom part) are not interchangeable. If one of the parts indicated becomes damaged then the complete decomposition vessel has to be sent to the factory for testing.

Accessories

Venting handle		C 7010.8
Quartz crucible		C 4
Ignition wire	(5 pieces)	C 5012.3
Cotton threads	(500 pieces)	C 710.4
Gelatin capsules	(100 pieces)	C 9
Acetobutyrate capsules	(100 pieces)	C 10
Benzoic acid (NBS)	(30 g)	C 43
Benzoic acid	(100 g)	C 43A
Benzoic acid tablets	(50 pieces)	C 723
Combustion bags 40x35	(100 pieces)	C 12
Combustion bags 70x40	(100 pieces)	C 12A
Oneway crucible		C 14
Foundation		C 5010.4

Other accessories can also be found in our **IKA-Catalogue**. The adress is included on the back of these operating instructions.


Use of oneway crucible

For the decomposition vessels AOD 1.1 are oneway crucibles C14 ease and shorten the preparation of combustion. They consist of a transparent plastik which burns completely to CO_2 and H_2O .

Advantages:

- The sample can be directly weighted in the oneway crucible.
- The weight of the oneway crucible C14 can be delivered from the weighing-machine to the calorimeter for the caculate the foreign energy.
- There's no more need for an cotton thread.
- The one way crucible C14 is a burning help. Samples burning incomplete in normal crucibles are burned better with the one way crucible C14. Thgere by better measuring results are reached.

Replacing of the crucible holder


For the use of the oneway crucible C14 the decomposition vessels AOD 1.1 must be put in the foundation C 5010.4.

This foundation C 5010.4 replaces the normal crucible holding.

Take thecover part and until the ignition wire on the crucible holder like discribed before.

Screw down the crucible holder from the shaft and replace it with the foundation C 5010.4.

Installation of the oneway crucible C14

The oneway crucible C14 must be put into the crucible holding of the foundation C 5010.4 so that the ignition wire touches the oneway crucible C14. It would the best to stick it under the ignition wire.

Technical data AOD 1.1


Operating pressure	max. 195 bar
Test pressure	280 bar
Volume	210 ml
Nuclear temperatures	$\geq 1300\text{ }^{\circ}\text{C}$
Operation temperatures	max. 50°C
Sample decomposition	isochor
Oxygen valve	non-returnable valve
(Vacuum operation not possible)	
Crucible holder	as in DIN 51 900
Material in contact with the product	
Lower part	2.4610.40
Cover	2.4610.40
Threated joint	1.4439

Guarantee

In accordance with IKA Terms and Conditions of Sale, the warranty period is 12 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

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IKA India Private Limited
IKA England LTD.
IKA Vietnam Company Limited
IKA Works (Thailand) Co. Ltd.