IKA UTL 25 ULTRA-TURRAX-Inline

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
● Scheduled Use

The ULTRA-TURRAX-INLINE UTL 25 is an inline dispersing instrument that operates according to the rotor-stator-principle. Its modular construction consists of flow chamber DK 25.11, ULTRA-TURRAX-drive T 25 S 5; S 6, dispersing element S 25 KV-25 G-IL and mounting support AD 25. With the help of this instrument, chemical and biotechnological processings up to + 220°C can be carried out.

The ULTRA-TURRAX-INLINE UTL 25 is suitable for inline-processings in which quickly reacting reaction partners are directly injected into the shearing head. It is also possible to directly immerse the chamber into liquids, thus making airfree dispersing with simultaneous circulating possible. In addition, the instrument can be used as a centrifugal pump.

● Unpacking

The ULTRA-TURRAX-INLINE UTL 25 is supplied in individual components. Please unpack these individual components carefully and look out for any damage. It is important to recognize a transportation damage already when unpacking. If this is the case, an immediate factual report is required (mail, rail or forwarding agency).

● Information on Application

With the help of the ULTRA-TURRAX-INLINE UTL 25, flowable and/or liquid media, emulsions and dispersions can be produced in the laboratory in continuous operation. Measurements of the dwelling time in the rotor/stator are just as possible as an operation under the exclusion of air.

The following set-ups are possible:


- continuous flow through
- circulation with pressure balance
- dosing of several liquids
- circulation of larger quantities

● Setting-up the Instrument

The dispersing instrument system ULTRA-TURRAX-INLINE UTL 25 is ready for use after fitting the individual components and connecting the threaded nozzle (inlet) and/or the threaded nozzle (outlet) to the reservoirs.

Assembly of the dispersing instrument system ULTRA-TURRAX-INLINE UTL 25:

- Loosen the 4 hexagonal screws on the flow chamber DK 25.11 and remove cover
- Insert dispersing element S 25 KV-25-IL
- Replace cover, check proper seat of O-ring, and tighten the 4 hexagonal screws
- Loosen ball grip screw 2 on mounting support AD 25
- Insert the unit flow chamber DK 25.11 and dispersing element into the mounting support AD 25
- Tighten ball grip screw 2
- Loosen ball grip screw 1 on flange of drive
- Push drive on dispersing shaft right to the stop, the on-off-switch is visible on the drive from the front
- Tighten ball grip screw 1 on flange of drive

Fitting the connecting hoses:

Threaded nozzle (inlet) LW 13 and threaded nozzle (outlet) LW 9

- pull hose clamp over hose
- push hose with hose clamp over threaded nozzle
- tighten hose clamp

● **Starting-Up**

Before starting-up the instrument, please check the conformity of the voltage indicated on the ID-plate of the drive with your mains voltage. If they do not agree, the instrument must not be operated.

 Separate operating instructions dispersing element S 25 KV-25 G-IL and ULTRA-TURRAX®-drive T 25 S 5 are attached for the individual components.

SAFETY WARNING:

The ball grip screws on the flange of the drive unit and on the mounting support might get loose by the vibration. For the sake of safety check the firm seat of the ball grip screws and retighten them, if necessary.

PLEASE NOTE!

Never operate the dispersing elements dry, as sealing and bearings will be destroyed if they are not cooled by the medium. The dispersing instrument system is not self-sucking. So please take care that the suction hose and the flow chamber are always filled with medium.


● Description of the Individual Components

Mounting Support AD 25

The mounting support is a sturdy light-metal cast construction. The mounting eye for the dispersing element is arranged in such a way that the entire ULTRA-TURRAX-INLINE UTL 25 unit stands safely on the laboratory bench, also without drive. The feet on the bottom side ensure a firm stand, even at maximum speed of the drive. With the help of the ball grip screw the dispersing element can be clamped to the mounting support.

Flow Chamber DK 25.11

The flow chamber incorporates the dispersing element and is provided with internal O-ring-sealings. The chamber casing houses the threaded nozzle LW 9 (outlet) and a hydraulic boring for the slide-ring-sealing. When supplied, the hydraulic boring is closed with a locking screw. When using the mode of operation „Continuous“-flow, it is recommended to establish a bypass-line into the feed hose with the help of the threaded nozzle M 5 that is included in the supply schedule. Any material not being processed is thus returned from the slide-ring-sealing area into the circuit.

Drive ULTRA-TURRAX T 25 S 5; S 6

Please see separate manual.

Dispersing element S 25 KV-25 G-IL

Please see separate manual.

● Maintenance, Cleaning and Sterilisation

Cleaning

Immediately after the dispersing instrument system has been operated, the dispersing element and the flow chamber must be cleaned so that the threads of rotor and stator do not stick together with remains of the substance. For this purpose the dispersing instrument system is operated in a solvent that removes the remains of the substance, yet does not attack the stainless steel. Because of the high velocity of flow, rotor and stator are already largely cleaned this way. For a thorough cleaning, however, the dispersing instrument system has to be dismantled.

Disassembly of the Dispersing Instrument System:

- Remove the hose connections
- Empty the flow chamber
- Loosen the ball grip screw 1 on flange of drive
- Draw off the drive from dispersing element
- Loosen hexagonal screws on cover of flow chamber and remove the cover
- Loosen ball grip screw 2 on mounting support, push flow chamber together with dispersing element out of mounting support
- Draw flow chamber carefully off from dispersing shaft
- Remove O-rings from cover and dispersing chamber

For Dismantling the Dispersing Element Please See Separate Manual!

Sterilisation

a) Chemical processing

By bactericidal solvents (formalin, phenol, alcohol) many cases of disinfection can be solved. It is important that the disinfectants are subsequently removed with sterilised water

b) Germ destruction by moist heat

This is carried out by saturated steam with an excess pressure of 120°C

c) Germ destruction by hot air

The hot air sterilisation is generally carried out at 160°C and/or 190°C.

The drive unit is maintenance-free, yet not wear-resistant (the motor carbons wear off during the course of time). In the interest of smooth operation and long durability the drive should be regularly cleaned.

For cleaning the surfaces of the drive and of the mounting support AD 25, exclusively

- a) water with a tenside containing detergent
- or b) isopropyl alcohol (only if very soiled) should be used.

SAFETY WARNING:

The drive may only be opened by a qualified technician – also in case of repair. Before opening it, the power supply plug has to be disconnected. Life parts inside the instrument can still be alive also a longer time after the plug has been disconnected.

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When ordering spare parts, please quote type of instrument, serial number as well as ordering number of the spare part.

Please submit only instruments for repair that have been cleaned and are free from substances dangerous to health.

● Materials Contacting the Product

Flow chamber	1.4301
Cover	1.4301
Threaded nozzles	1.4301
Locking screw	1.4301
O-rings	FFKM

● Technical Data

Quantity of flow	lt/min	5 ... 15
max temp. of medium	°C	120
threaded nozzle inlet		LW 13
internal-Ø of inlet-hose	mm	13
threaded nozzle outlet		LW 9
internal-Ø of outlet-hose	mm	9

● Accessories

S 25 KV-25 G-IL	dispersing element
DK 25.11	flow chamber
AD 25	mounting support
DZM 5	IKA-TRON revolution counter
RLS 8 x 33	measuring sensor
R 1822	plate stand
R 1823	plate stand
R 181	boss head
DTM 11	digital temperature measuring
PT 100.3	instrument
R 340	temperature measuring sensor,
	stainless steel
	IKA-safety-reservoir

● Guarantee

You have acquired an Original IKA-Laboratory Instrument which meets the highest requirements of quality and design. In accordance with the IKA-terms of warranty, this instrument is guaranteed for a period of 12 months. In the event of any claim under guarantee, please contact your laboratory dealer.

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