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Safety instructions

- **Read the assembly instructions completely before starting up and follow the safety instructions.** If this accessory is used with another device, observe also its operating instructions.
- Refer also to the safety instructions of operating instructions of the rotary evaporator.
- Keep the assembly instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations. **When working under a vacuum in particular!**
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. Otherwise there is a risk of:
 - splashing liquids,
 - body parts, hair, clothing and jewellery getting caught,
 - injury as a result of glass breakage.
- Ensure that there is sufficient space above the device as the glass assembly may exceed the height of the device.
- Prior to each use, always check the device, accessories and especially the glass parts for damage. Do not use damaged components.
- Ensure that the glass assembly is tension-free! Danger of cracking as a result of:
 - stress due to incorrect assembly,
 - external mechanical hazards,
 - local temperature peaks.
- Operation with excess pressure is not permitted (cooling water pressure 0.5 bar, minimum cooling water flow rate 30 l/h, maximum cooling water flow rate 100 l/h).
- The glass equipment must always be ventilated when working under normal pressure (e.g. open outlet at cooler) in order to prevent a pressure build-up.
- Evacuated glass vessels must not be heated on one side; the evaporating flask must rotate during the heating phase.
- The glassware is designed for operation under a vacuum of up to 1 mbar. The equipment must be evacuated prior to heating (see operating instructions of the rotary evaporator, chapter "Commissioning"). The equipment must only be aired again after cooling. When carrying out vacuum distillation, undensified vapours must be condensed out or safely dissipated. If there is a risk that the distillation residue could disintegrate in the presence of oxygen, only inert gas must be admitted for stress relief.
- **Danger of burning!** The heating bath, tempering medium, evaporator piston and glass assembly can become hot during operation and remain so for a long time afterwards! Let the components cool off before continuing work with the device.
- A vacuum may form inside the glassware in case of power outage of the rotation drive.
- Imbalance may result in uncontrolled resonance behaviour of the rotation drive or the assembly. Glass apparatus may be damaged or destroyed. In the event of unbalance or unusual noises, switch off the rotation drive immediately or reduce the speed.
- **Caution!** The surface of the evaporating flask is slippery when silicone oil is used as a tempering medium in the heating bath!

Product information

The glassware is made of borosilicate glass 3.3.
The coating of the glassware coated is made out of the material SURLYN.
The glassware coated by SURLYN are suitable for the use in a temperature range from -30 upto +60 °C.


Delivery scope

Item	Designation	Quantity glassware diagonal	Quantity glassware vertical
1	Receiving flask	1	1
2	Clamp NS 29, stainless steel	1	1
4	Stopcock	1	1
5	Tube	-	1
6	Condenser	1 Diagonal condenser	1 Vertical condenser
7	Connection	1 Introduction sleeve	1 Vacuum connection
8	Clamp NS 29, plastic	1	1
10	Evaporator piston 1000 ml	1	1
11	Ball joint clamp RV 05.10	1	1
12	Screw joint cap	4	4
13	Hose connection	4	4
18	Washer	-	1

Mounting vacuum seal (illustrations zoomed in see manual rotary evaporator)

Fit the condenser seal

- ☞ Place the condenser seal in the condenser receptacle and fit the glassware to the device according to the assembly instructions.

First use


- ☞ Insert vapor tube (1).
- ☞ Insert the vacuum seal (2). Pay attention to the correct positon of the vacuum seal!
- ☞ Slide the union nut (3b) over the flange on the condenser (3a).
- ☞ Also slide the annular spring (3c) over the flange on the condenser (3a).

- ☞ Position the condenser (3a) on top of the seal (2).


- ☞ Tighten the union nut (3b) by hand onto the threaded flange. Run the device for 20 minutes at 120 rpm. Finally, retighten the union nut (3b) by hand.

- ☞ Secure the vertical glassware with the Velcro® band.

Note: See operating instructions for the rotary evaporator!

Mounting the washer (illustrations zoomed in see manual rotary evaporator)

Note:
Pay attention to the correct positon of the washer.
The PTFE tube (5) with the washer (18) can be mounted optionally on the vertical glass sets.
It serves for filling the evaporator flask when there is a vacuum in the glass set.
Solvent can be sucked through the PTFE tube into the evaporator flask by opening the stopcock (4).

Cleaning

The seal on the glass condenser should be checked at regular intervals and replaced if necessary.

For cleaning disconnect the mains plug!

To clean the glassware use only water with a detergent that contains tensides, or use isopropylalcohol for stubborn soiling.

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