


Femoral

Interlocking System



● Femoral Nail Features

- Stainless steel 316 LVM.
- Can be used as a minimal invasive technique by cannulated instruments or by traditional technique.
- Universal use - left or right.
- Cannulated.
- Anatomical Nail.
- With Transverse proximal locking screws.
- With dynamic proximal hole.
- With accurate distal target device for distal screws.
- Self tapping locking screws.
- Complete instrument set with the newly designed instruments like
 - Cannulated awl
 - Smart reduction device (Finger like)
 - Newly designed gripper to help the surgeon in manipulation with the guide wire
 - Mechanical Distal target device instrument set



Reconstruction nail



Retrograde nail

EgiFix medical provides several solutions for the femoral fixation with intramedullary nails that can meet the surgeons need according to the type of fracture and the patients' need

Reconstruction nail

with screws to catch the neck of the femur
different entry port through the greater trochanter
two bowing to meet the anatomical bowing of the femur


Retrograde nail

can be an advantage when used for distal femur fracture type A and C according to AO classification

Femoral Interlocking Instrument Set



Awl Cannulated(EGTI2-000-01)



Awl Curved(EGTI2-000-02)




Gripper(EGTI2-000-03)



Flexible Reamers(EGTI2-000-04)



Proximal Assembly(EGTI2-000-05)



Distal target measure(EGTI2-000-06)



Distal target Device 1(EGTI2-000-07)



Stabilizing Rod 6.0mm(EGTI2-000-08)



Metal Ruler(EGTI2-000-09)

Drill Guide 4.0mm(EGTI2-000-10)
Drill Guide 5.2mm(EGTI2-000-11)

Drill Bit 4.0mm(EGTI2-000-12)



Reamer 6.0mm(EGTI2-000-13)



Depth Gauge(EGTI2-000-14)



Reduction Device(EGTI2-000-15)



Sliding Hummer(EGTI2-000-16)



Extractor(EGTI2-000-17)



Angled Trocar(EGTI2-000-18)

Screw Driver(EGTI2-000-19)
Trocar 8.0mm(EGTI2-000-20)

8 Extension wally El Ahd St, Hadaek Elkobba - Cairo - Egypt.

Tel.: +(202) 2601 9519


E-mail:info@egifix.com

website:www.egifix.com

Femoral Interlocking Instrument Set



Awl Cannulated(EGTI2-000-01)



Awl Curved(EGTI2-000-02)




Gripper(EGTI2-000-03)



Flexible Reamers(EGTI2-000-04)



Proximal Assembly(EGTI2-000-05)



Distal target measure(EGTI2-000-06)



Distal target Device 1(EGTI2-000-07)



Stabilizing Rod 6.0mm(EGTI2-000-08)



Metal Ruler(EGTI2-000-09)

Drill Guide 4.0mm(EGTI2-000-10)
Drill Guide 5.2mm(EGTI2-000-11)

Drill Bit 4.0mm(EGTI2-000-12)



Reamer 6.0mm(EGTI2-000-13)



Depth Gauge(EGTI2-000-14)



Reduction Device(EGTI2-000-15)



Sliding Hummer(EGTI2-000-16)



Extractor(EGTI2-000-17)



Angled Trocar(EGTI2-000-18)

Screw Driver(EGTI2-000-19)
Trocar 8.0mm(EGTI2-000-20)

8 Extension wally El Ahd St, Hadaek Elkobba - Cairo - Egypt.

Tel.: +(202) 2601 9519

E-mail:info@egifix.com

website:www.egifix.com