


# LEAD PLACEMENT - STEP 1

- **V<sub>1</sub> - R Sternal  
Border 4th ICS**
- **V<sub>2</sub> - L Sternal  
Border 4th ICS**
- **V<sub>4</sub> - L Mid-  
clavicular 5th ICS**





# LEAD PLACEMENT - STEP 2

- **V<sub>3</sub>** - Between V<sub>2</sub> & V<sub>4</sub>
- **V<sub>6</sub>** - L Mid-Axillary at Same Horizontal Level as V<sub>4</sub>
- **V<sub>5</sub>** - L Anterior-Axillary at same horizontal level as V<sub>4</sub> & V<sub>6</sub>

