


IKA

designed for scientists

IKA BC 1000



Technische Information
Technical information


DE
EN

Allgemeine Informationen

Der Becherglasdeckel **IKA BC 1000** kann mit 1000 ml- und 600 ml-Bechergläsern zum Dispergieren und Rühren verwendet werden. Der Becherglasdeckel verhindert ein Spritzen und senkt den Geräuschpegel beim Rühren und Dispergieren.

Der Deckel **BC 1000** ist aus Silikon und gegen die meisten der bei Versuchen benutzten Säuren und Basen beständig.

Installation



Hinweis: Während der Anwendung muss die Lippe des Becherglasdeckels am Ausguss des Becherglases dicht anliegen.

❶ bei 1000 ml-Becherglas.

❷ bei 600 ml-Becherglas.

Eigenschaften und Vorteile



Material

Silikon



Max. Temperatur

200 °C



Spülmaschinenfest

ja



Vorteil

kein Spritzen



Vorteil




senkt den Geräuschpegel

General information

IKA BC 1000 beaker cap can be used with 1000 ml and 600 ml beakers for dispersing and stirring experiments. The beaker cap can prevent splashing, as well as reduce the noise level while stirring and dispersing.

The BC 1000 beaker cap is made of silicone. It can resist most acidic and alkaline media used in experiments.

Installation



Note: The lip on the beaker cap must fit the beaker mouth tightly during application.

❶ for 1000 ml beaker.

❷ for 600 ml beaker.

Character and advantage

	material
	silicone

	max. temperature
	200 °C

	dishwasher safe
	yes

	advantage
	no splashing

	advantage
	reduce noise



designed for scientists

IKA-Werke GmbH & Co. KG

Janke & Kunkel-Straße 10, 79219 Staufen, Germany

Phone: +49 7633 831-0, Fax: +49 7633 831-98

eMail: sales@ika.de

USA

IKA Works, Inc.

Phone: +1 910 452-7059

eMail: usa@ika.net

KOREA

IKA Korea Ltd.

Phone: +82 2 2136 6800

eMail: sales-lab@ika.kr

BRAZIL

IKA Brazil

Phone: +55 19 3772 9600

eMail: sales@ika.net.br

MALAYSIA

IKA Works (Asia) Sdn Bhd

Phone: +60 3 6099-5666

eMail: sales.lab@ika.my

CHINA

IKA Works Guangzhou

Phone: +86 20 8222 6771

eMail: info@ika.cn

POLAND

IKA Poland Sp. z o.o.

Phone: +48 22 201 99 79

eMail: sales.poland@ika.com

JAPAN

IKA Japan K.K.

Phone: +81 6 6730 6781

eMail: info_japan@ika.ne.jp

INDIA

IKA India Private Limited

Phone: +91 80 26253 900

eMail: info@ika.in

UNITED KINGDOM

IKA England LTD.

Phone: +44 1865 986 162

eMail: sales.english@ika.com

VIETNAM

IKA Vietnam Company Limited

Phone: +84 28 38202142

eMail: sales.lab-vietnam@ika.com

Discover and order the fascinating products of IKA online:

www.ika.com



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide

Technical specifications may be changed without prior notice.

20000004065_25000611b_BC 1000_092020