

PRODUCT INSTRUCTIONS

(For LD777-1)


WARNING!

- This chair has a 250 lbs weight capacity.
- Hang cushion over chair frame to dry if it gets wet to prevent mold/mildew.

Remove all protective packing material and position base on floor. Adjust all the foot levelers level, and put on the cushion, attached the cushion on the chair frame. Now the chair is ready for using!

Warning: DO NOT use acid or alkaline washing liquid to clean the product.

A		CHAIR	1
B		FOOT LEVELER	4
C		SEAT CUSHION	1



This parts have been assemble by factory already!